

BLICKMAN



Bassinets

Transport infants safely and comfortably

Includes:


- Bassinet Kits
- Bassinet Frames
- Basket & Mattress



THE BLICKMAN DIFFERENCE

- Our bassinets are renowned for their reliable 304 stainless steel, welded construction, durability and optimal performance which has made them the preferred choice among maternity nurses for decades.
- Our inverted “U-frame” bassinet design provides a stable and comfortable resting place for newborns. Our standard 4” ball bearing swivel casters, with two brakes, ensure safe transport.
- We offer various storage choices such as shelves, drawers and cabinets. Cabinets can be closed for more security or open for easier access.
- Our clear acrylic infant basket includes a cardholder for baby identification and rests in specially formed support rods so it can be set in the flat or Trendelenburg position.

Key Product Features



Bassinets

Bassinet Kits

We offer bassinet kits so you can get everything you need in one order. Our bassinet kits include:

- ✓ The stainless steel bassinet frame of your choosing
- ✓ The acrylic basket
- ✓ The foam mattress pad



Now Included

BASSINET KITS

All of our bassinets are the same dimensions:

Width: 31" x Height: 41 1/2" x Depth: 17 3/4"



Description:

- Upper shelf has four sided guardrail to protect supplies
- Two shelves are 8" and 23" from the floor
- Shelf dimensions (WxHxD): 31" x 17 3/4"

Description:

- Storage drawer that is counter weighted for stability
- One shelf, 8" from the floor
- Drawer inside dimensions (WxHxD): 19 3/4" x 6 1/2" x 13 7/8"

Description:

- Two 17" deep storage compartments with doors
- Right side has removable shelf
- Cabinet inside dimensions (WxHxD): 29" x 13 3/4" x 15 7/8"

Description:

- Two open 17" deep storage compartments
- Right side includes non-removable shelf
- Cabinet inside dimensions (WxHxD): 29" x 13 3/4" x 16 5/8"

Name	Bassinet Kit with Open Shelves	Bassinet Kit with Drawer	Bassinet Kit with Closed Cabinet	Bassinet Kit with Open Cabinet
Model #	8044SS-KIT	8047SS-KIT	8048SS-KIT	8049SS-KIT
Item #	1118044KIT	1118047KIT	1118048KIT	1118049KIT
Shipping Weight	69 lbs.	98 lbs.	98 lbs.	77 lbs.

Custom Designs Available Upon Request
Please email us at info@blickman.com

Bassinets

Parts Ordered Separately

BASSINETS

All of our Bassinets are the same dimensions:

Width: 31" x Height: 38 1/2" x Depth: 17 3/4"



Name	Bassinet with Open Shelves	Bassinet with Drawer	Bassinet with Closed Cabinet	Bassinet with Open Cabinet
Model #	8044SS	8047SS	8048SS	8049SS
Item #	1118044000	1118047000	1118048000	1118049000
Ship Weight	53 lbs.	83 lbs.	82 lbs.	61 lbs.

ACCESSORIES



Basket

Dimensions: 31 1/8" x 9 7/8" x 16 1/8"

Seamless, clear acrylic bassinet basket with cardholder for easy identification



Mattress

Dimensions: 26" x 2" x 12 3/4"

2" thick foam padding with staph check cover. Flame retardant

Name	Basket	Mattress
Model #	3070P	3069MP
Item #	9841307001	8441306901
Ship Weight	12 lbs.	3 lbs.

Custom Designs Available Upon Request

Please email us at info@blickman.com

Bassinets

Construction Specifications

Construction Specifications

Universal Dimensions **Bassinets:** Width: 31" x Height: 38 1/2" x Depth: 17 3/4"
Bassinet Kits: Width: 31" x Height: 41 1/2" x Depth: 17 3/4"

Materials 304 stainless steel

Gauge Legs: 16GA
Shelves: 18 GA
Body Parts and Removable Inside Shelves: 20GA

Welding All welding by heliarc method and exposed welds are ground smooth

Finish Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish

Casters 4" ball bearing, swivel casters, two with brakes

Notes

- All models utilize specially formed support rods so basket can be set in the flat or Trendelenburg position
- Infant basket and mattress pad can be sold separately

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Custom Designs Available Upon Request

Please email us at info@blickman.com

Toll Free: (800) 247-5070
Local: (973) 330-0557

BLICKMAN

www.blickman.com
info@blickman.com

1099 Wall Street West, Suite 250A, Lyndhurst, NJ 07071