

BLICKMAN



Exam/ Patient Room

Improve Workflow and Efficiency

Includes:

- Folding Utility Cart
- Overbed Table
- Mayo Stands
- Foot Stools
- Exam Stools
- Patient Chairs



THE BLICKMAN DIFFERENCE

- Exam rooms are typically the first place in which patients seek medical care, except in the event of an emergency.
- We provide a wide variety of exam room furniture products to fit your needs. From overbed tables to folding utility carts to vinyl exam stools and patient chairs.
- Our stainless steel equipment provide sterile and long lasting solutions that will withstand many patients rotating through your exam room.
- If you're looking for more storage solutions, go to www.blickman.com or view our cabinet catalog.

Key Products



Folding Utility Cart



Overbed Table



Mayo Stands



Exam Stools



Patient Chairs



Foot Stools



Exam Room

Utility Cart & Overbed Table

FOLDING UTILITY CART

Description: Our folding utility cart is versatile and can assist with almost any application such as movement of mail, patient records, and laboratory materials. Choose any combination of 6" or 12" deep wire baskets to complete your cart. The metal frame is designed for easy assembly and folds for storage. Rubber stoppers support and protect the wire racks to provide a quiet ride on four, 3" hard rubber casters.

Folding Utility Cart Dimensions:

Width: 29 1/4" x Height: 38 1/2" x Depth: 19 1/2"



Name	Model #	Item #	Shipping Weight
Folding Utility Cart	2440	2422440000	21 lbs.
ORDER SEPARATELY			
6" Wire Basket	2441	2422441000	5 lbs.
12" Wire Basket	2442	2422442000	6 lbs.

OVERBED TABLE

Description: Our overbed table comes in an attractive mark resistant laminate top in a walnut wood grain. It moves up or down with a simple one touch adjustment. Security features include a built-in column lock and four point base. The table travels smoothly on four 2" hooded casters.

Overbed Table Dimensions

Width: 30" x Height: 30" (max) x Depth: 15 1/2"



Name	Overbed Table
Model #	3400
Item #	2833400000
Shipping Weight	27 lbs.

Exam Room

Mayo Stands

CHROME MAYO STANDS

Description: Each of our chrome Mayo stands contain a 12 5/8" x 19 1/8" stainless steel tray and are height adjustable using a Tru-Loc™ friction knob. We offer different base and pole set ups to meet your varying needs.

Middle Pole

Base: Five leg cast aluminum base with wall savers™ and 2" twin wheel hooded casters



Double Pole

Base: Chrome loop base with two 2" twin wheel hooded back casters



Single Pole

Base: Two pronged chrome base with single hard rubber back casters



Name	Middle Pole	Double Pole	Single Pole
Height Adjustment	29 1/2" - 49"	34 1/2" - 54"	34 1/2" - 54"
Model #	1501	1510	1515
Item #	0661501000	0661510000	0661515000
Shipping Weight	14 lbs.	18 lbs.	12 lbs.

Exam Room

Mayo Stands

NEWARK MAYO STANDS

Description: Each of our stainless steel Mayo stands contain a 12 5/8" x 19 1/8" stainless steel tray and are height adjustable using a Tru-Loc™ friction knob. We offer different base set ups to meet your varying needs.

Lightweight

Base: Two pronged stainless steel base with two hard rubber back casters



Standard



Base: Two pronged stainless steel base with two roller blade back casters



Four Caster

Base: Two pronged stainless steel base with four roller blade casters



Name	Lightweight	 Standard	 Four Caster
Height Adjustment	37 1/4" - 64"	39 1/2" - 62"	39 1/2" - 62"
Model #	1520SS	8841SS	8841MR
Item #	0671520000	0638841000	0638841200
Shipping Weight	26 lbs.	36 lbs.	36 lbs.

Exam Room




Mayo Stands

BENJAMIN MAYO STAND

Description: Our Benjamin Mayo stand design is made for ease and durability. Our state of the art double post design minimizes lateral movement and increases overall stability. We offer this design in three different tray sizes.



Via a foot operated mechanism, you can adjust the height from 39 1/2" to 62" off the ground. A guard over the foot prevents accidental release.

Name	Benjamin Mayo Stands	
Tray Size	12 5/8" x 19 1/8" (Standard)	
Model #	8867SS	
Item #	0668867000	
Shipping Weight	40 lbs.	
Tray Size	16" x 21"	
Model #	8869SS	
Item #	0668869000	
Shipping Weight	41 lbs.	
Tray Size	20" x 25"	
Model #	8871SS	
Item #	0668871000	
Shipping Weight	43 lbs.	

Exam Room

Foot Stools

8" HIGH FOOT STOOLS

Description: Our MR conditional approved Kent foot stools include one step 8" off the ground. They are industry standard, fully welded, stainless steel construction with non-slip rubber treads and feet tips. 18" and 30" width sizes are available. Depth is 12". Include a 41" tall handrail for additional support.

8" stools have the same dimensions:
Width: 18" or 30" x Height: 8" x Depth: 12"



**18" Wide
Kent Stool**

**30" Wide
Kent Stool**

Name	18"W	30"W
With Handrail (41" tall)		
Model #	7757MR-HR	7761MR-HR
Item #	1017757100	1017761100
Shipping Weight	20 lbs.*	30 lbs.*
Without Handrail		
Model #	7757MR	7761MR
Item #	1017757000	1017761000
Shipping Weight	11 lbs.	20 lbs.

*Shipping weight with an asterisk is the combined shipping weight of separate parts

Exam Room

Foot Stools

9" HIGH FOOT STOOLS

Description: Our 9" high foot stools come in either chrome or stainless steel and consist of one step 9" off the ground with non-slip rubber treads and feet tips. A 35 3/4" tall handrail for additional support is available for the chrome stool.

All 9" stools have the same dimensions:

Width: 14 1/4" x Height: 9" x Depth: 11 3/8"



Stainless Steel Stool

Name	Model #	Item #	Ship Weight
Stainless Steel Stool	1260SS	1011260000	8 lbs.

Parts Also Sold Separately



Chrome Stool with Handrail



#1250 - Stool



#651 - Handrail

Name	Model #	Item #	Ship Weight
Chrome Stool with HR	1251	1011251000	11 lbs.*

ORDER SEPARATELY

Chrome Stool	1250	1021251640	7 lbs.
Handrail (35 3/4" tall)	651	7900651000	4 lbs.

*Shipping weight with an asterisk is the combined shipping weight of separate parts

Exam Room

Foot Stools

TWO STEP FOOT STOOLS

Description: Our MR conditional approved two step foot stools are industry standard, fully welded, stainless steel construction with H-braces for additional reinforcement, non-slip rubber treads for secure footing, and feet tips for firm floor placement.



Step Surface Areas
 Top Step: 36" x 12"
 Bottom Step: 36" x 24"

Heights of Steps
 Top Step: 15 3/4"
 Bottom Step: 8"

Name	2-Step Foot Stool
With Handrails (41" high)	
Model #	7764MR-HR
Item #	1012041096
Shipping Weight	65 lbs.

Two Step Foot Stool with Side Handrails



Step Surface Areas
 Top Step: 24" x 8"
 Bottom Step: 24" x 16"

Heights of Steps
 Top Step: 15 3/4"
 Bottom Step: 8"

Name	The Donnelly
With Handrails (41" high)	
Model #	7763MR-HR
Item #	1017763100
Shipping Weight	40 lbs.*
Without Handrail	
Model #	7763MR
Item #	1017763000
Shipping Weight	27 lbs.

The Donnelly with Handrail

Exam Room

Foot Stools

OUR EZ Stacking Stool allows you to:

- ✓ Create wider platforms by locking side-to-side
- ✓ Create higher platforms by stacking top-to-bottom
- ✓ Store efficiently and transport easily with our corresponding dolly

EZ STACKING STOOLS





Description: Our EZ stacking stools and dolly are made with our durable stainless steel. The stools are covered with a non-slip rubber tread and feet tips.



EZ Stool Dimensions:
20" x 6 1/4" x 15"



Dolly Dimensions:
24 3/4" x 39" x 17 5/16"

Name	 EZ Stacking Stool		Dolly
Model #	8862SS	8862MR	8862SS-D
Item #	1018862000	10188620MR	1018862300
Ship Weight	14 lbs.	14 lbs.	20 lbs.

Exam Room

Exam Stools

4 LEG SPINDLE EXAM STOOLS

Description: We offer two four-leg spindle adjusted exam stools with the option of 2" casters or rubber feet tips. Choose our Passaic stool and sit comfortably on a 15" diameter, 3" thick polyfoam padded seat upholstered in vinyl or choose our Clifton stool for an easy-to-clean stainless steel seat with a non-slip vinyl insert. Add a removable pad slipcover to your Clifton stool for a more comfortable option.



Spindle

Adjusts 12"
Rubber Tips: 19" - 31"
Casters: 21" - 33"



Name

CLIFTON

PASSAIC

Standard

With Pad

Standard

Casters

Model #	7745SS-C	7745SS-PC	7714-PSSC
Item #	1027745001	1027745003	1027714001
Shipping Weight	17 lbs.	17 lbs.	17 lbs.

Feet

Model #	7745SS	7745SS-P	7714-PSS
Item #	1027745000	1027745002	1027714000
Shipping Weight	16 lbs.	16 lbs.	16 lbs.



Options: Available in standard black. Consult factory for requests in other colors for large order quantities

Exam Room

Exam Stools

5 LEG SPINDLE EXAM STOOL

Description: Our five-leg spindle exam stool rests on a 24" diameter base. Choose between a black composite base or an aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.



Spindle

Adjusts 5 7/8"
Height: 20 7/8" - 26 3/4"



Name	5 Leg Spindle Stool	
	Aluminum Base	5 Leg Spindle Stool with backrest
Model #	1113	1113WB
Item #	1021113000	1021113200
Shipping Weight	22 lbs.	26 lbs.
	Black Composite Base	
Model #	1113BLK	1113WB-BLK
Item #	1021113100	1021113300
Shipping Weight	22 lbs.	26 lbs.

Options: Foot ring, 75 different vinyl colors - see upholstery catalog

Exam Room

Exam Stools

FOOT OPERATED EXAM STOOL

Description: Our five-leg foot operated exam stool have a push ring under a 24" diameter aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.



Foot Pedal

Adjusts 7"
Height: 18 3/8" - 25 3/8"



Name	Foot Operated Stool	Foot Operated Stool with Backrest
	Aluminum Base	
Model #	1205	1206
Item #	1041205000	1041206000
Shipping Weight	22 lbs.	28 lbs.

Options: Foot ring, 75 different vinyl colors - see upholstery catalog

Exam Room

Exam Stools

HAND OPERATED EXAM STOOLS

Description: Our five-leg hand operated exam stool rests on a 24" diameter base. Choose between a black composite base or an aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.

Standard Stools



Hand Lift

Adjusts 7"
Height: 18 3/8" - 25 3/8"



Name	Hand Operated Stool	Hand Operated Stool with Backrest
	Aluminum Base	
Model #	1210	1210WB
Item #	1041210000	1041210200
Shipping Weight	22 lbs.	26 lbs.
	Black Composite Base	
Model #	1210BLK	1210WB/BLK
Item #	1041210100	1041210300
Shipping Weight	22 lbs.	26 lbs.

Options: Foot ring, 75 different vinyl colors - see upholstery catalog



Exam Room

Exam Stools

HAND OPERATED EXAM STOOLS

Economy Stools

Description: Our economy hand operated exam stool rests on a 24" diameter black composite base with black plastic twin wheel casters. Sit in comfort on a 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. this model is available with and without a back.

Adjusts 7" from 18 3/8" to 25 3/8" off the ground



Name	Economy Stool	
	Without Back	With Back
Model #	1200	1201
Item #	1041200125	1041201325
Shipping Weight	22 lbs.	26 lbs.

Options: Foot ring, 75 different vinyl colors - see upholstery catalog

Deluxe Stools

Description: Our deluxe hand operated exam stool comes standard with a foot ring for leg-joint comfort. These stools rest on a 24" diameter black composite or aluminum base, which comes with Wall Saver™ tips to protect against scuff marks and offer a textured surface for foot grip. Both bases travel easily on 2" hard rubber casters. Sit in comfort on a 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding with a backrest.

Adjusts 7" from 21 3/8" to 28 3/8" off the ground.



Name	Deluxe Stool	
	Aluminum Base	Black Composite Base
Model #	1212	1212BLK
Item #	1041212025	1041212125
Shipping Weight	29 lbs.	29 lbs.

Options: 75 different vinyl colors - see upholstery catalog

Foot Ring as an Accessory: Compatible on all 5 leg stools

Part #	Description	Shipping Weight
7901243000	Foot Ring for Leg Joint Comfort	6 lbs.



Exam Room

Chairs

PATIENT CHAIRS

Description: Our patient chairs are designed for comfortable seating in patient exam rooms. Our standard size is 24" x 33" x 22 5/8" and come with or without arm rests. All our patient chairs are available in 75 vinyl colors.



Name	Chair	Chair with Arms
Model #	1130	1130WA
Item #	1051130025	1051130125
Shipping Weight	21 lbs.	17.5 lbs.

Options: 75 vinyl colors. See upholstery catalog or visit www.blickman.com/chairs for all available colors

Toll Free: (800) 247-5070
Local: (973) 330-0557



www.blickman.com
info@blickman.com

1099 Wall Street West, Suite 250A, Lyndhurst, NJ 07071